


Nasal Spray Pump

Multi-dose Nasal Spray Pump



Optional Nasal Actuator



Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410 for screw-on
20mm for snap-on
20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging

Fully physical & chemical validation, high level microbiological safety

Comply with CP, EP, USP food & drug direct contact regulations

US DMF# 031363

Oral Spray Pump

Multi-dose Oral Spray Pump



Optional Oral Actuator

L



O



J



I



SR



LR



3"



N



Sf1



Es1



Sf2



LF



Op1



Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410

for screw-on

20mm for snap-on

20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging

Fully physical & chemical validation, high level microbiological safety

Comply with CP, EP, USP food & drug direct contact regulations

US DMF# 031364

Topical spray Pump

Multi-dose Topical Spray Pump

Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410 for screw-on
20mm for snap-on
20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging

Fully physical & chemical validation, high level microbiological safety

Comply with CP, EP, USP food & drug direct contact regulations

Optional Topical Actuator



Ear Spray pump

Multi-dose Ear Spray Pumps

Optional Ear Actuator

Sf1



Es1



Sf2



LF



Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410 for screw-on

20mm for snap-on

20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging

Fully physical & chemical validation, high level microbiological safety

Comply with CP, EP, USP food & drug direct contact regulations



Closures

Screw on



18/410



18/415



18Din



18GL



18TE

Dosage (mcl)	45	50	65	67	70	78	85	90	100	110	120	130	140	150
	170	180	200	220										



20/410



20/415

Dosage (mcl)	45	50	65	67	70	78	85	90	100	110	120	130	140	150
	170	180	200											



24/410

Dosage (mcl)	100	143	170	
-----------------	-----	-----	-----	--

Closures

Snap on



15mm

Dosage
(mcl)

80



20mm

Dosage
(mcl)

45
200

50

65

67

70

78

85

90

100

110

120

130

140

150

170

180

Crimp on



18mm



20mm/G1



20mm/G2

Dosage
(mcl)

45
200

50

65

67

70

78

85

90

100

110

120

130

140

150

170

180

Innovation & Patented Products



Child Resistant Nasal Pump

Key Features:

Comply with 16CFR1700.20

Tamper-evident closure

15ml and 30ml HDPE bottle are available

Invention Patented in China, USA, Europe and Korea



Preservative Free Nasal Pump

Key Features:

Air filter membrane

Air filter compatible with medicine solution

0.14ml dosage available

Invention Patented in China, USA, Europe and Korea



Preservative Free Eye Drop Bottle

Key Features:

Air filter membrane

Passed microbiological challenging test

5ml, 7ml, 10ml volume available